



Eradikit[®] BVDV 1 and 2 discrimination test



Features

- Indirect ELISA based on recombinant E2 antigens
- Serological discrimination between BVDV 1 and 2 pestivirus
- Estimation of herd protection vs BVDV1 and 2 following vaccination
- Spot test vs BVDV1 and 2
- Applicable to serum, plasma and individual or bulk milk
- 5 plates (230samples)



Product name	Code	Plates	Samples
Eradikit [®] BVDV1-2	2838-5	5F8	230

Information:

www.in3diagnostic.com

info@in3diagnostic.com

Order:

order@in3diagnostic.com